

METHOD CONTROL PARAMETERS

Method Information For: C:\MSDCHEM\1\METHODS\TFMPP.M
Method Sections To Run:

- () Save Copy of Method With Data
- () Instrument Control Pre-Run Cmd/Macro =
- () Data Analysis Pre-Run Cmd/Macro =
- (X) Data Acquisition
- (X) Data Analysis
- () Instrument Control Post-Run Cmd/Macro =
- () Data Analysis Post-Run Cmd/Macro =

Method Comments:

END OF METHOD CONTROL PARAMETERS

Base Conversion → 2 methods

- ① Add sodium Hydroxide to partial tablet in tube.
Soak overnight. Then extract with hexane

or

- ② Add drop of conc Ammonium Hydroxide to partial
tablet in a vial (tube) and then add petroleum ether

* HP-1ms gives best results (@ 25m x 0.2mm x 0.33mm)

* HP-5ms gives okay separation (~~0-30~~ 30m x 0.25mm x 0.33mm)

ASD

11-17-10

C:\MSDCHEM\1\METHODS\TFMPP.M
Wed Nov 17 11:34:14 2010

Control Information

Sample Inlet : GC
Injection Source : GC ALS
Mass Spectrometer : Enabled

=====

6890 GC METHOD

OVEN

Initial temp: 120 'C (On) Maximum temp: 325 'C
Initial time: 2.50 min Equilibration time: 0.25 min
Ramps:
Rate Final temp Final time
1 10.00 280 0.00
2 0.0 (Off)
Post temp: 0 'C
Post time: 0.00 min
Run time: 18.50 min

FRONT INLET (SPLIT/SPLITLESS)

BACK INLET (UNKNOWN)

Mode: Split
Initial temp: 250 'C (On)
Pressure: 11.61 psi (On)
Split ratio: 40:1
Split flow: 39.9 mL/min
Total flow: 43.5 mL/min
Gas saver: On
Saver flow: 20.0 mL/min
Saver time: 2.00 min
Gas type: Helium

COLUMN 1

COLUMN 2

Capillary Column
Model Number: Agilent 19091S-433
HP-5MS 5% Phenyl Methyl Siloxane
Max temperature: 325 'C
Nominal length: 30.0 m
Nominal diameter: 250.00 um
Nominal film thickness: 0.25 um
Mode: constant flow
Initial flow: 1.0 mL/min
Nominal init pressure: 11.61 psi
Average velocity: 38 cm/sec
Inlet: Front Inlet
Outlet: MSD
Outlet pressure: vacuum

(not installed)

FRONT DETECTOR (NO DET)

BACK DETECTOR (NO DET)

SIGNAL 1

SIGNAL 2

Data rate: 20 Hz
Type: test plot
Save Data: Off
Zero: 0.0 (Off)

Data rate: 20 Hz
Type: test plot
Save Data: Off
Zero: 0.0 (Off)

Range: 0
Fast Peaks: Off
Attenuation: 0

Range: 0
Fast Peaks: Off
Attenuation: 0

COLUMN COMP 1
(No Detectors Installed)

COLUMN COMP 2
(No Detectors Installed)

THERMAL AUX 2

Use: MSD Transfer Line Heater

Description:

Initial temp: 295 'C (On)

Initial time: 0.00 min

#	Rate	Final temp	Final time
1	0.0(Off)		

POST RUN

Post Time: 0.00 min

TIME TABLE

Time	Specifier	Parameter & Setpoint
------	-----------	----------------------

GC Injector

Front Injector:

Sample Washes	3
Sample Pumps	3
Injection Volume	1.00 microliters
Syringe Size	10.0 microliters
PreInj Solvent A Washes	1
PreInj Solvent B Washes	1
PostInj Solvent A Washes	1
PostInj Solvent B Washes	1
Viscosity Delay	0 seconds
Plunger Speed	Fast
PreInjection Dwell	0.00 minutes
PostInjection Dwell	0.00 minutes

Back Injector:

No parameters specified

Column 1 Inventory Number : AB001

Column 2 Inventory Number :

MS ACQUISITION PARAMETERS

General Information

Tune File : stune.u
Acquisition Mode : Scan

MS Information

Solvent Delay : 1.50 min
EM Absolute : False
EM Offset : 0
Resulting EM Voltage : 1423.5

[Scan Parameters]

Low Mass : 25.0
High Mass : 500.0
Threshold : 200
Sample # : 2 A/D Samples 4
Plot 2 low mass : 40.0
Plot 2 high mass : 500.0

[MSZones]

MS Source : 230 C maximum 250 C
MS Quad : 150 C maximum 200 C

END OF MS ACQUISITION PARAMETERS

TUNE PARAMETERS for SN: US65125812

Trace Ion Detection is OFF.

EMISSION : 34.610
ENERGY : 69.922
REPELLER : 19.904
IONFOCUS : 64.533
ENTRANCE_LE : 0.000
EMVOLTS : 1423.529
AMUGAIN : 1646.000
AMUOFFSET : 129.438
FILAMENT : 2.000
DCPOLARITY : 0.000
ENTLENSOFFS : 14.557@ 3 14.557@ 50 9.788@ 69 10.541@131 10.541@219 13.051@
414 13.302@502 13.302@1049
MASSGAIN : -797.000
MASSOFFSET : -39.000

END OF TUNE PARAMETERS

END OF INSTRUMENT CONTROL PARAMETERS

DATA ANALYSIS PARAMETERS

Method Name: C:\MSDCHEM\1\METHODS\TFMPP.M

Percent Report Settings

Sort By: Retention Time

Output Destination

Screen: No
Printer: Yes
File: No

Integration Events: Meth Default

Generate Report During Run Method: No

Signal Correlation Window: 0.020

Qualitative Report Settings

Peak Location of Unknown: Apex minus Start of Peak

Library to Search Minimum Quality

C:\Database\SLI.L 80
C:\Database\NIST05a.L 75
C:\Database\PMW_TOX2.L

Integration Events: Meth Default

Report Type: Summary

Output Destination

Screen: No
Printer: Yes
File: No

Generate Report During Run Method: No

Quantitative Report Settings

Report Type: Summary

Output Destination

Screen: Yes
Printer: No
File: No

Generate Report During Run Method: No

Calibration Last Updated: Mon Dec 17 08:44:05 2007

TFMPP.M Wed Nov 17 11:34:02 2010

Reference Window: 10.00 Percent
Non-Reference Window: 5.00 Percent
Correlation Window: 0.02 minutes
Default Multiplier: 1.00
Default Sample Concentration: 0.00

Compound Information

1) mass 284 ()

Ret. Time 7.280 min., Extract & Integrate from 6.780 to 7.780 min.

Signal	Rel Resp.	Pct. Unc.(rel)	Integration
Tgt 283.90			*** METH DEFAULT ***

Lvl ID	Conc ()	Response
3	1.000	-1
4	10.000	-1
5	100.000	-1
2	0.100	-1

Qualifier Peak Analysis ON
Curve Fit: Avg. RF

2) mass 283 ()

Ret. Time 7.280 min., Extract & Integrate from 6.780 to 7.780 min.

Signal	Rel Resp.	Pct. Unc.(rel)	Integration
Tgt 283.90			*** METH DEFAULT ***

Lvl ID	Conc ()	Response
3	1.000	-1
4	10.000	-1
5	100.000	-1
2	0.100	-1

Qualifier Peak Analysis ON
Curve Fit: Avg. RF

3) mass 247 ()

Ret. Time 7.280 min., Extract & Integrate from 6.780 to 7.780 min.

Signal	Rel Resp.	Pct. Unc.(rel)	Integration
Tgt 246.90			*** METH DEFAULT ***

Lvl ID	Conc ()	Response
3	1.000	-1
4	10.000	-1
5	100.000	-1
2	0.100	-1

Qualifier Peak Analysis ON
Curve Fit: Avg. RF

4) mass 212 ()

Ret. Time 7.280 min., Extract & Integrate from 6.780 to 7.780 min.

Signal	Rel Resp.	Pct. Unc.(rel)	Integration
Tgt 212.00			*** METH DEFAULT ***

Lvl	ID	Conc ()	Response
3		1.000	-1
4		10.000	-1
5		100.000	-1
2		0.100	-1

Qualifier Peak Analysis ON
Curve Fit: Avg. RF

END OF DATA ANALYSIS PARAMETERS

Wed Nov 17 11:34:26 2010